

OLSC MS10a (Magnetic Stirrer 10 Liters)

When looking for magnetic stirrers, users look for a device which delivers reliable performance and years of durability, but most manufacturers tend to oversell and under deliver which results in a frustrating experience for the user. Magnetic stirrers come with promises of exceptional performance, quality components and high durability, but in just a short time, the reality becomes apparent. Users face grave issues like the top getting corroded or the motor is not powerful enough to stir large quantities of solutions, its speed remains low etc.

With detailed insights about the requirements and needs of users and their feedback of existing stirrers, the R&D department at Neuation came up with the design of OLSC MS10a magnetic stirrer to develop a product that delivers high performance without compromising on flexibility or durability. With its Plastic Integrated SS Top surface, the device is extremely resistant to corrosion while its long lasting DC motor runs at top speed without needing much of a maintenance. The single knob operation makes the machine extremely convenient to use. The same motor is capable enough to stir large quantities of solutions at once for higher productivity and less wastage of time and effort.



Product Features



DC Motor For Quiet Performance And A Long Life.



Stirring Capacity Up To 10 Litres



Variable Speed Setting From 400-2200 RPM.



Plug & play



Small Footprint Saves Valuable Bench Space And Also Makes It Suitable For Use In Hoods And Cold Room.

No corrosive Plastic Integrated SS Top increases life of instrument



Non-slip Strength Rubber Feet For Maximum Stability Even At High-speeds.



IP42 Class Compliant Rating For Unmatched Durability.



Safe Low Voltage Device



Specifications

Specifications	OLSC-MS10a
Capacity	10 Litres
Motor Type	DC
Variable Speed	400 - 2200 RPM
Set-up plate dimensions	135 mm
Weight	1.5 kg
Dimension (WxDxH)	248 x 156 x 92 mm
Protection class	IP 42
Set-up plate material	Plastic Integrated SS Top Plate

Ordering Information

OLSC-MS10a		
Product Name	Cat #	Details
OLSC-MS10a	HS20-2614 EU/US/AU/UK	Magnetic Stirrer

Standard Accessories	
Part Name	Cat #
Stir Bar	B1

