


## STUDENT COMPOUND MICROSCOPE



### SALIENT FEATURE

Traditional Horse shoe base, with body inclinable at 90°, with mechanical tube length of 160mm.

Double lens bright field ABBE Condenser.

Swing-in and swing-out filter holder

Separate coarse and fine focusing mechanism

Eyepiece Huygenian eyepiece 10x & 15x

Nosepiece: Triple revolving Nosepiece

Objectives JIS 10X, 40X/ 45X (S), 100X (Oil, S)

Mechanical Stage Stage plate size 120 x 120mm with clips for holding the slide

Illumination Illumination system Plano concave reflector at the base (fork mounted).