

Mini Vortex Mixer



Stepless speed regulation
Wide speed range



High quality and colorful
silicone



Press to start / Compact
body



Brushless motor / Durable

TECHNICAL DATA

Model	DSVM-300S
Control Mode	Knob
Motor	Brushless motor
Amplitude	4.5mm (Circular oscillation)
Maximum Sample Processing Capacity	50ml
Speed Range	0-3200rpm
Shell Protection Level	IP42
Power	12W
Voltage	100-240V, 50/60Hz
Product Size	100*93*72mm
Allowable Ambient Temperature	5-40°C,80%RH
Gross / Net weight	0.56/0.37kg