



Plasma Expressor Xpress



Xpress is specially designed plasma expressor for separating out the blood components. The front panel of the

Xpress is loaded with high tension spring for quick transfer of plasma into other bag.

Xpress is ergonomically designed for long lasting trouble free operations

Salient Features

- User friendly with see through front panel to allow visibility of layer movement
- High tension spring loaded thick compression plate for fast plasma expression
- Ergonomic handle with auto adjustable lever for easy operation

Technical Specifications

Parameter	Unit	OLSC-Xpress
Dimensions (W x D x H)	mm	165x230x280
Weight	kg	3.5
Body Construct		Powder Coated CRCA Steel
Front Panel		Acrylic