



AH

Automatic 5-Diff Hematology Analyzer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Automatic 5-Diff Hematology Analyzer

DSK-8550

5-Diff Differential of WBC
Double Chamber of WBC and RBC/PLT
29 parameters + 4 Histogram + 1
Scattergram
Up to 60 Sample/hour



- Compatible with two methods of automatic injection and manual open injection. The automatic mode supports uninterrupted and continuous addition, and emergency samples can be inserted at any time. First detection
- Fully automatic sample injection completely simulates manual sample inversion and mixing.
- Intelligent WBC floating landmark technology, floating classification technology can automatically adapt to the changes of scatter plots, effectively identify the boundaries of various types of cells, and obtain more accurate white blood cell classification results.
- Automatic detection of blocked holes and blockage removal.

Model	DSK-8550
Sample capacity	50 supports emergency analysis
Function	360 rotating bar cod scanning
Printer	Can be connected to an external laser printer
Display	14 inches color-full touch screen display
Storage	Auto save more then 200.000 analysis result analysis
Data transpiration	USBX4, Rs232, LAN, HDMI supports multilingual settings
Dimension	70X64X67 cm
Weight	57 kg